

Perfusion Neuroimaging:

**Cerebral Blood Flow
through Brain Tissue Matter**


Recap: MR Angiography-MRA

BLOOD FLOW VELOCITY in Artery

Time-of-Flight(TOF) MRA: “Bright Blood Imaging”

- Magnetic resonance angiography is a group of techniques based on magnetic resonance imaging to image blood vessel.
- MRA is often used to evaluate the arteries of the neck and brain, the thoracic and abdominal aorta, the renal arteries, and the legs.
- Time of Flight(TOF) MRA was the dominant non-contrast bright-blood method for imaging the human vascular system.



Time-of-flight MRA at the level of Circle of Willis



Recap: MRA.

MRI Velocitometry:

Color Coding of Blood Flow Velocity in cm/s across an arterial cross-section




Time Variation with Arterial Pulsing:

Velocity Profile in Carotid Artery at different time points as Heart Expands & Contracts

PERFUSION

- Perfusion is the rate of blood flow through capillary circulation of an organ or tissue.
- Measurement of perfusion requires the use of tracer molecules or particles that may be intravascular, extracellular or freely diffusible.

Measurement of Cerebral Blood Flow with **PET or MRI** (Arterial Spin Labeling - ASL)




- Uses magnetically labeled arterial blood water as an endogenous flow tracer
- Potentially provide quantifiable CBF in classical units (mL/min per 100 gm of tissue)

TYPES OF **PERFUSION** MRI

- There are 3 types of perfusion using MRI.


1) Dynamic Susceptibility Contrast

Most widely used ,rapid **T2*** imaging after gadolinium bolus.




2) Dynamic Contrast Enhanced

T1-imaging after gadolinium bolus; measure vascular permeability.



3) Arterial Spin Labeling-



Uses magnetically tagged water in blood (not Gadolinium) as tracer.




ARTERIAL SPIN LABELING

- Arterial spin labeling (ASL) is an Magnetic Resonance Imaging method for measuring perfusion using the patient's own water molecules as tracers.
- ASL is a non-invasive , no exogenous contrast injection required.
- ASL is reliable and reproducible, repeated measurements possible.
- Absolute quantification for blood flow.
- **Note:**
Unit of Perfusion or Cerebral Blood Flow is say **X ml. of Blood** flowing **through 100 gram** of Brain Tissue **per minute.**

Recap: 180° Inversion Pulse: DTI b-value




ARTERIAL SPIN LABELING PRINCIPLES




ARTERIAL SPIN LABELING PRINCIPLES

- First “control ” images are acquired through area of interest (cerebral hemispheres).
- Next, “tagged” pulse are applied to a slab of tissue proximally from the imaging volume that **inverts the magnetization** of water molecules in this slab
- In next couple of seconds, most of these “magnetically labeled” molecules lying within the vessels will flow into imaging volume.
- **These tagged water molecules exchange their magnetization with those in the static tissue.** The area of interest is re-imaged
- Tagged images are subtracted from control image and the **result of this subtraction is perfusion-weighted image.**



PARAMETERS OF ASL

- Blood flow(BF), Blood Volume(BV) and Mean transit time are basic parameters.



- Blood flow and Blood volume are assumed to be constant.

BASIL TOOL BOX FOR PERFUSION MAPPING



Basil

Bayesian Inference for Arterial Spin Labelling

Input Data | Structure | Calibration | Distortion Correction | Analysis

Data contents

Input Image (None) 

Number of PLDs 1 - +

Number of repeats <Unknown>

Data order

Grouping order Label/Control pairs ▾

☒ Label/Control pairs Label then control ▾

Acquisition parameters

Labelling cASL/pcASL ▾


Bolus duration (s) Constant  1.80 - +

Bolus durations (s) 1.8

PLDs (s) 1.8

Readout 3D (eg GRASE) ▾ Time per slice (ms)  10.00 - +

Data preview - perfusion weighted image



Use scroll wheel to change slice, double click to change view

Data order preview



Basil

Bayesian Inference for Arterial Spin Labelling

Input Data | Structure | Calibration | Distortion Correction | **Analysis**

Data preview - perfusion weighted image

Basic analysis options

Output Directory

Browse

☐ Brain Mask

(None)



☐ Analysis which conforms to 'White Paper' (Alsop et al 2014)

Initial parameter values

Arterial Transit Time (s)

1.30

-

+

T1 (s)



1.30

-

+

T1b (s)

1.65

-

+

Inversion Efficiency

0.85

-

+

Analysis Options

☒ Adaptive spatial regularization on perfusion

☐ Incorporate T1 value uncertainty

☐ Include macro vascular component

☒ Fix label duration


☐ Partial Volume Correction

Use scroll wheel to change slice, double click to change view

Data order preview


Illustrative Example - 1

GRAPHICAL ANALYSIS OF CEREBRAL BLOOD FLOW AND DIFFERENT AGE GROUP



Analysis of CBF with ageing (Normal Subjects):
There is no statically significant change

Illustrative Example - 2



Analysis of CBF between **Normal subjects and Alzheimer's patients** of the Same Age Groups:

There is statistically significant decrease.

CONCLUSION OF THE STUDY

- From observations we have concluded that value of cerebral blood flow in normal patients is more than the Alzheimer disease patients.